

# DCAF12 Rabbit pAb

货号: B17683

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Rat                                                                                            |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 50kDa                                                                                                |
| 实测分子量 | 47kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HepG2,Rat kidney                                                                                     |
| 细胞定位  | Cytoplasm,centrosome,cytoskeleton,microtubule organizing center                                      |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human DCAF12 (NP_056212.1). |
| 序列   | MARKVVSRRKAPASPGAGSDAQGPQFGWDHSLHKRRLPPVKRSLVYYLKNREVRLQNETSYSRVLHGYYAAQQLPSLL KEREFH                             |

靶点信息

|      |                                                                                                                                                                                                                                     |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a WD repeat-containing protein that interacts with the COP9 signalosome, a macromolecular complex that interacts with cullin-RING E3 ligases and regulates their activity by hydrolyzing cullin-Nedd8 conjugates. |
| 基因ID | 25853                                                                                                                                                                                                                               |

|       |                                    |
|-------|------------------------------------|
| 基因名   | DCAF12                             |
| Swiss | Q5T6F0                             |
| 别名    | DCAF12;CT102;KIAA1892;TCC52;WDR40A |

产品验证



Western blot analysis of DCAF12 expressed in HepG2,Rat kidney using DCAF12 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)