

GPR82 Rabbit pAb

货号: B17686

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	38kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Mouse heart,Rat heart
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPR82 (NP_543007.1).
序列	MNNNTTCIQPSMISSMALPIIYILLCIVGVFGNTLSQWIFLTKIGKKTSTHIYLSHLVTANLLVCSAMPFMSIYFLKGFQWEYQS AQCRVVNFLGTLSMH

靶点信息

研究背景	The protein encoded by this gene is an orphan G protein-coupled receptor of unknown function. The encoded protein is a member of a family of proteins that contain seven transmembrane domains and transduce extracellular signals through heterotrimeric G proteins.
基因ID	27197

基因名	GPR82
Swiss	Q96P67
别名	GPR82

产品验证



Western blot analysis of GPR82 expressed in Mouse testis,Mouse heart,Rat heart using GPR82 Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug p er lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s

实验步骤

访问官网浏览详情: www.ablybio.cn