

CHDH Rabbit pAb

货号: B17709

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	65kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HepG2,Mouse thymus,Mouse kidney,Mouse liver,Mouse pancreas,Mouse testis,Rat kidney
细胞定位	Mitochondrion inner membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 245-594 of human CHD H (NP_060867.2).
序列	LHPALSRTNLKAEATLVSRLVFEGTRAVGVEYVKNQSHRAYASKEVILSGGAINSPQLMLSGIGNADDLKKLGIPVVCH LPGVGQNLQDHLIYIQQACTRPITLHSAQKPLRKVCIGLEWLWKFTGEGATAHLETGGFIRSQPGVPHPDIQFHFLPSQVI DHGRVPTQQEAYQVHVGPMTGTSVGWLKLRSANPDHPVIQPNYLSTETDIEDFRLCVKLTREIFAQEALAPFRGKELQP GSHIQSDKEIDAFVRAKADSAYHPSCTCKMGQPSDPTAVVDPQTRVLGVENLRVVDASIMPSMVSGNLNAPTIMIAEKAA DIIKGQPALWDKDVVPYKPRTLATQR

靶点信息

研究背景	The protein encoded by this gene is a choline dehydrogenase that localizes to the mitochondrion. Variations in this gene can affect susceptibility to choline deficiency. A few transcript variants have been found for this gene, but the full-length nature of only one has been characterized to date.
------	---

基因ID	55349
基因名	CHDH
Swiss	Q8NE62
别名	CHDH

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn