

HES6 Rabbit pAb

货号: B17711

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	12kDa/23kDa/24kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human HES6 (NP_061115.2).
序列	MAPPAAPGRDRVGREDEDGWETRGRKARKPLVEKKRRRARINESLQELRLLLAGAEVQAKLENAEVLTLVRRVQGVLRGRAREREQLQAEASERFAAGYIQCMHEVHTFVSTCQ AIDATVAAELLNHLLSMP LREGSSFQDLLGDALAGPPRAPGRSGWPAGGAPGSPGPSPPGPGDDLCS DLEEAPEAELSQA PAEGPDLVPAALGSLTTAQIARSVWRPW

靶点信息

研究背景	This gene encodes a member of a subfamily of basic helix-loop-helix transcription repressors that have homology to the Drosophila enhancer of split genes. Members of this gene family regulate cell differentiation in numerous cell types. The protein encoded by this gene functions as a cofactor, interacting with other transcription factors through a tetrapeptide domain in its C-terminus. Alternatively spliced transcript variants encoding different isoforms have been described.
------	---

基因ID	55502
基因名	HES6
Swiss	Q96HZ4
别名	HES6;C-HAIRY1;HES-6;bHLHb41;bHLHc23

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn