

GPR108 Rabbit pAb

货号: B17716

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:5000
理论分子量	60kDa
实测分子量	76KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,U2OS,Mouse colon,Rat thymus
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 33-260 of human GPR108 (NP_001073921.1).
序列	RIHQLALTGEKRADIQLNSFGFYTNLSLEVELSVLRLGLREAEKSLLVGFSLSRVRSGRVRSYSTRDFQDCPLQKNSSSFLVFLINTKDLQVVRKYGEQKTLFIFPGLLPEAPSKPGLPKPQATVPRKVDGGGTSAASKPKSTPAVIQGPSGKDKDLVLGLSHLNNSYNFSFHVIGSQAEEGQYSLNFHNCNNSVPGKEHPFDITVMIREKNPDGFLSAAEMPL

靶点信息

研究背景	Predicted to be involved in negative regulation of toll-like receptor signaling pathway. Located in trans-Golgi network.
基因ID	56927

基因名	GPR108
Swiss	Q9NPR9
别名	GPR108;LUSTR2

产品验证



Western blot analysis of GPR108 expressed in Mouse colon,Rat thymus,U-87MG using GPR108 Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug p er lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

实验步骤

访问官网浏览详情: www.ablybio.cn