

LNx2 Rabbit pAb

货号: B17757

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	76kDa
实测分子量	76kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-375,A-549,Rat thymus
细胞定位	plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human LNx2 (NP_699202.1).
序列	MGTTSDEMVSVEQTSSSSSLNPLCFECGQQHWTRENHLYNQNEVDDDLVCHICLQPLLQPLDTPCGHTFCYKCLRNFLQ EKDFCPLDRKRLHFKLCKKSSILVHKLLDKLLVLCPFSSVCKDVMQRCdleAHLKNRCPGASHRRVALERRKTSRTQAEIEN ENGPTLLDPAGTlsPEADCLGTGAVPVERHLTSASLSTWSEEPGLDNPafEESAGADTTQQPLSLPEGEITTIeIHRSN

靶点信息

研究背景	Predicted to enable ubiquitin-protein transferase activity. Predicted to be involved in protein ubiquitination.
基因ID	222484

基因名	LNX2
Swiss	Q8N448
别名	LNX2;PDZRN1

产品验证



Western blot analysis of LNX2 expressed in A-375,A-549,Rat thymus using LNX2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn