

Cpsf4l Rabbit pAb

货号: B17764

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	20kDa
实测分子量	22kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Mouse brain
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of mouse Cpsf4l (NP_080958.2).
序列	MEEVIAGLQGFTFAFELDVESQKGTGLLPFQGMKSSSAVCNFFAKGLCVKGMLCPLRHEQGEKLVVCKHWLRGLCRKS DCCDFLHQYDVSKMPVCYFHSKFGNCSNKECLFLHLKPVLKLQDCPWYNQGFCKEVGPLCKYRHVHQVLCPNYFTGFCP EGPQCQFGHPKMSPPFHPSNVKLQPVNQPWDSTTSTGNTLVSPFQAKPMVHGQRRWSLPQACSSRAHPAP

靶点信息

研究背景	Predicted to be involved in pre-mRNA cleavage required for polyadenylation. Predicted to be part of mRN A cleavage and polyadenylation specificity factor complex. Is expressed in central nervous system and re tina. Orthologous to human CPSF4L (cleavage and polyadenylation specific factor 4 like).
基因ID	52670

基因名	
Swiss	A6NMK7
别名	CPSF4L;0610010C04Rik;1500000C01Rik;D11Ert636e;Cpsf4l

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn