

MPP7 Rabbit pAb

货号: B17768

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	42kDa/65kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung
细胞定位	Cell junction,Membrane,Peripheral membrane protein,adherens junction,tight junction
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 297-576 of human MPP7 (NP_775767.2).
序列	RLALRRPEILVQPLKVSNRKSSGFRKSFRLSRKDKKTNKSMYECKKSDQYDTADVPTYEEVTPYRRQTNEKYRLVVLGVPV GVGLNELKRKLLISDTQHYGVTPHTTRARRSQESDGYEYIFISKHLFETDVQNNKFIEYGEYKNNYYGTSIDSVRSVLAKN KVCLLDVQPHTVKHLRTLLEFKPYVIFIKPPSIERLRETRKNAKIISSRDDQGAAPFTEEDFQEMIKSAQIMESQYGHLFDKIII NDDLTVAFNELKTTFDKLETETHWVPVSWLHS

靶点信息

研究背景	The protein encoded by this gene is a member of the p55 Stardust family of membrane-associated guanylate kinase (MAGUK) proteins, which function in the establishment of epithelial cell polarity. This family member forms a complex with the polarity protein DLG1 (discs, large homolog 1) and facilitates epithelial cell polarity and tight junction formation. Polymorphisms in this gene are associated with variations in site-specific bone mineral density (BMD). Alternative splicing results in multiple transcript variants.
------	---

基因ID	143098
基因名	MPP7
Swiss	Q5T2T1
别名	MPP7

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn