

Thromboxane synthase Rabbit pAb

货号: B17894

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200
理论分子量	52kDa/60kDa/65kDa
实测分子量	55KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse lung,Mouse spleen,Rat lung
细胞定位	Endoplasmic reticulum membrane,Multi-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 407-533 of human Thromboxane synthase (NP_001052.3).
序列	AFRFTREAAQDCEVLGQRIPAGAVLEMAVGALHHDPEHWSPETFNPERFTAEARQQHRPFTYLPFGAGPRSCLGVRLGL LEVKLTLLHVLHKFRFQACPETQVPLQLESKSALGPKNGVYIKIVSR

靶点信息

研究背景	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. However, this protein is considered a member of the cytochrome P450 superfamily on the basis of sequence similarity rather than functional similarity. This endoplasmic reticulum membrane protein catalyzes the conversion of prostglandin H2 to thromboxane A2, a potent vasoconstrictor and inducer of platelet aggregation. The enzyme plays a role in several pathophysiological processes including hemostasis, cardiovascular disease, and stroke. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	6916
基因名	TBXAS1
Swiss	P24557
别名	TBXAS1;BDPLT14;CYP5;CYP5A1;GHOSAL;THAS;TS;TXAS;TXS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn