

PTPRA Rabbit pAb

货号: B17906

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	89kDa/90kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,22Rv1,BT-474,Mouse brain,Mouse kidney
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 20-142 of human PTPRA (NP_543031.1).
序列	NNATTVAPSVGITRLINSSTAEPVKEEAKTSNPTSSLTSLSVAPTFSPNITLGPTYLTTVNSSDSDNGTTRTASTNSIGITISPN GTWLDPNQFTDARTEPWEGNSSTAATTPETFPPSDETP

靶点信息

研究背景	The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains an extracellular domain, a single transmembrane segment and two tandem intracytoplasmic catalytic domains, and thus represents a receptor-type PTP. This PTP has been shown to dephosphorylate and activate Src family tyrosine kinase s, and is implicated in the regulation of integrin signaling, cell adhesion and proliferation. Three alternativ ely spliced variants of this gene, which encode two distinct isoforms, have been reported.
基因ID	5786
基因名	PTPRA
Swiss	P18433
别名	PTPRA;HEPTP;HLPR;HPTPA;HPTPalpha;LRP;PTPA;PTPRL2;R-PTP-alpha;RPTPA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn