

PRDM1 Rabbit pAb

货号: B17909

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	76kDa/87kDa/91kDa
实测分子量	100KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 70-210 of human PRDM 1 (NP_001189.2).
序列	GADGGTSVQAEASLPRNLLFKYATNSEEVIGVMSKEYIPKGTRFGPLIGEITYNDTVPKNANRKYFWRIYSRGELHHFIDGF NEEKSNWMRYVNPAAHSPREQNLAAQCNGMNIYFYTIKPIPANQELLVWYCRDFAERLHY

靶点信息

研究背景	This gene encodes a protein that acts as a repressor of beta-interferon gene expression. The protein bind s specifically to the PRDI (positive regulatory domain I element) of the beta-IFN gene promoter. Transcript ion of this gene increases upon virus induction. Two alternatively spliced transcript variants that encode d ifferent isoforms have been reported.
基因ID	639

基因名	PRDM1
Swiss	O75626
别名	PRDM1;BLIMP1;PRDI-BF1

产品验证



Western blot analysis of PRDM1 expressed in Jurkat,Mouse liver using PRDM1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn