

DLC1 Rabbit pAb

货号: B17927

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	50kDa/54kDa/114kDa/122kDa/126kDa/170kDa
实测分子量	171kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung
细胞定位	Cell junction,Cytoplasm,Membrane,Peripheral membrane protein,focal adhesion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 204-463 of human DLC1 (NP_079043.3).
序列	DAPKVNAVDTLNVKDIAPEKQLLNSAVIAQRRKPDPPKDENERSTCNVVQNEFLDTPCTNRGLPLLKTDFGSCLLQPPSC PNGMSAENGLEKSGFSQHQNKSPPKVKAEDGMQCLQLKETLATQEPTDNQVRLRKRKEIREDRDRARLDSMVLLIMKLD QLDQDIENALSTSSPSGTPTNLRRHVPDLESGSESGADTISVNQTRVNLSSDTESTDLPSSTPVANSGTKPKTTAIQGISE KEKAGKLTFWFCFLANLF

靶点信息

研究背景	This gene encodes a GTPase-activating protein (GAP) that is a member of the rhoGAP family of proteins w hich play a role in the regulation of small GTP-binding proteins. GAP family proteins participate in signalin g pathways that regulate cell processes involved in cytoskeletal changes. This gene functions as a tumor suppressor gene in a number of common cancers, including prostate, lung, colorectal, and breast cancers . Multiple transcript variants due to alternative promoters and alternative splicing have been found for thi s gene.
------	--

基因ID	10395
基因名	DLC1
Swiss	Q96QB1
别名	DLC1;ARHGAP7;HP;STARD12;p122-RhoGAP

产品验证



Western blot analysis of DLC1 expressed in Mouse lung using DLC1 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn