

PADI4 Rabbit pAb

货号: B17932

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:20 - 1:50
理论分子量	74kDa
实测分子量	74kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,SKOV3
细胞定位	Cytoplasm,Cytoplasmic granule,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 290-663 of human PADI 4 (NP_036519.2).
序列	VFRVAPWIMTPNTQPPQEVYACSFENEDFLKSVTTLAMKAKCKLTICPEEENMDDQWMQDEMEIGYIQAPHKTLPPVVFDS PRNRGLKEFPPIKRVMGPDFGYVTRGPQTGGISGLDSFGNLEVSPPVTVRGKEYPLGRILFGDSCYPSNDSRQMHQALQDF LSAQVQAPVKLYSDWLSVGHVDEFLSFVPAPDRKGFRLLASPRSCYKLFQEQQNEGHEALLFEGIKKKKQQKIKNLS NKTLRHNSFVERCIDWNRELLKRELGLAESDIIDIPQLFKLKEFSKAEAFFPNMVNMLVLGKHLGIPKPFPGPVINGRCCLEEK VCSLLEPLGLQCTFINDFFTYHIRHGEVHCGTNNVRKPFSEFKWWNMVP

靶点信息

研究背景	This gene is a member of a gene family which encodes enzymes responsible for the conversion of arginin e residues to citrulline residues. This gene may play a role in granulocyte and macrophage development l eading to inflammation and immune response.
------	---

基因ID	23569
基因名	PADI4
Swiss	Q9UM07
别名	PADI4;PAD;PAD4;PADI5;PDI4;PDI5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn