

LIPG Rabbit pAb

货号: B17937

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa/56kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,A-549,Mouse liver
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 251-500 of human LIPG (NP_006024.1).
序列	GCGLNDVLGSIAYGTITEVVKCEHERAVHLFVDSLVDKPSFAFQCTDSNRFKKGICLSRKNRCNSIGYNAKKMRNKR NSKMYLKTRAGMPFRVYHYQMKIHVFSYKNMGEIPTFYVTLYGTNADSQTLPLEIVERIEQNATNTFLVYTEEDLGDLLKIQ LTWEGASQSWYNLWKEFRSYLSQPRNPGRELNIRRIRVKSGETQRKLTFC TEDPENTSISPGRELWFRKCRDGWRMKNE TSPTVELP

靶点信息

研究背景	The protein encoded by this gene has substantial phospholipase activity and may be involved in lipoprote in metabolism and vascular biology. This protein is designated a member of the TG lipase family by its se quence and characteristic lid region which provides substrate specificity for enzymes of the TG lipase fam ily.
------	--

基因ID	9388
基因名	LIPG
Swiss	Q9Y5X9
别名	LIPG;EDL;EL;PRO719

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn