

MLLT4 Rabbit pAb

货号: B17951

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	
实测分子量	200-250kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Mouse lung
细胞定位	adherens junction,cell junction,cell-cell contact zone,cell-cell junction,cytosol,nuclear speck,nucleoplasm, plasma membrane,tight junction
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1510-1680 of human MLLT4 (NP_001353249.1).
序列	SLSPDPWKRAKEKLEKQQMHIVDMLSKEIQELQSKPDRSAEESDRLRKLMLEWQFQKRLQESKQKDEDEDEEEDDD VDTMLIMQRLEAERRARLQDEERRRQQLEEMRKREAEDRARQEEERRRQEEERTKRDAEEKRRQEEGYYSRLEAERRRQHDEAARRLLEPEA

靶点信息

研究背景	This gene encodes a multi-domain protein involved in signaling and organization of cell junctions during embryogenesis. It has also been identified as the fusion partner of acute lymphoblastic leukemia (ALL-1) gene, involved in acute myeloid leukemias with t(6;11)(q27;q23) translocation. Alternatively spliced transcript variants encoding different isoforms have been described for this gene, however, not all have been fully characterized.[provided by RefSeq, May 2011]
------	---

基因ID	4301
基因名	AFDN
Swiss	P55196
别名	AF6; MLL-AF6; MLLT4; I-afadin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn