

# IL10RA Rabbit pAb

货号: B17957

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:100 - 1:500 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	63kDa
实测分子量	90-110KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant protein of human IL10RA.
序列	

## 靶点信息

研究背景	The protein encoded by this gene is a receptor for interleukin 10. This protein is structurally related to int erferon receptors. It has been shown to mediate the immunosuppressive signal of interleukin 10, and thu s inhibits the synthesis of proinflammatory cytokines. This receptor is reported to promote survival of pro genitor myeloid cells through the insulin receptor substrate-2/PI 3-kinase/AKT pathway. Activation of this receptor leads to tyrosine phosphorylation of JAK1 and TYK2 kinases. Two transcript variants, one protein -coding and the other not protein-coding, have been found for this gene.
------	--

基因ID	3587
基因名	IL10RA
Swiss	Q13651
别名	IL10RA;CD210;CD210a;CDW210A;HIL-10R;IL-10R1;IL10R

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)