

ERCC3 Rabbit pAb

货号: B17995

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	89kDa
实测分子量	89kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Mouse testis,Mouse thymus,Mouse brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 513-782 of human ERCC3 (NP_000113.1).
序列	PMSPEFYREYVAIKTKKRILLYTMNPNKFRACQFLIKFHERRNDKIIVFADNVFALKEYAIRLNKPYYGPTSQGERMQILQNFK HNPKINTIFISKVGDTSFDLPEANVLIQISSHGGSSRRQEAQRLGRVLRAKKGMVAEEYNAFFYSLVSQDTQEMAYSTKRQR FLVDQGYSFKVITKLAGMEEEDLAFSTKEEQQLLKVLAATDLDAEEEVVAGEFGSRSSQASRRFGTMSSMSGADDTV YMEYHSSRSKAPSKHVVHPLFKRFRK

靶点信息

研究背景	This gene encodes an ATP-dependent DNA helicase that functions in nucleotide excision repair. The encoded protein is a subunit of basal transcription factor 2 (TFIIH) and, therefore, also functions in class II transcription. Mutations in this gene are associated with Xeroderma pigmentosum B, Cockayne's syndrome, and trichothiodystrophy. Alternative splicing results in multiple transcript variants.
------	--

基因ID	2071
基因名	ERCC3
Swiss	P19447
别名	ERCC3;BTF2;GTF2H;RAD25;TFIIH;TTD2;XPB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn