

PEDF/SERPINF1 Rabbit pAb

货号: B18005

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	46kDa
实测分子量	46KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse liver,Human serum
细胞定位	Melanosome,Secreted
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human PEDF/SERPINF1 (NP_002606.3).
序列	IKSSFVAPLEKSYGTRPRVLTGNPRLDLQEINNWWQAQMKGKLARSTKEIPDEISILLGVAHFKGQWVTKFDSRKTSLEDF YLDEERTVRVPMMSDPKAV

靶点信息

研究背景	The protein encoded by this gene is a member of the serpin family, although it does not display the serin e protease inhibitory activity shown by many of the other serpin family members. The encoded protein is secreted and strongly inhibits angiogenesis. In addition, this protein is a neurotrophic factor involved in n eural differentiation in retinoblastoma cells.
基因ID	5176

基因名	SERPINF1
Swiss	P36955
别名	SERPINF1;EPC-1;OI12;OI6;PEDF;PIG35

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn