

Sall4 Rabbit pAb

货号: B18028

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	65kDa/112kDa
实测分子量	75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,Mouse liver,Rat liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 853-1053 of human Sall 4 (NP_065169.1).
序列	PGMTPLLAQPRRQAKQHGCTRCGKNFSSASALQIHERTHTGEKPFVCNICGRAFTTKGNLKVHYMTHGANNNNSARRGR KLAIENTMALLGTDGKRVSEIFPKEILAPSVNVDPVVWNQYTSMLNGLAVKTNEISVIQSGGVPTLPVSLGATSVVNNATV SKMDGSQSGISADVEKPSATDGVPKHQFPHFLEENKIAVS

靶点信息

研究背景	This gene encodes a zinc finger transcription factor thought to play a role in the development of abducen s motor neurons. Defects in this gene are a cause of Duane-radial ray syndrome (DRRS). Alternative splici ng results in multiple transcript variants encoding different isoforms.
------	---

基因ID	57167
基因名	SALL4
Swiss	Q9UJQ4
别名	SALL4;DRRS;HSAL4;ZNF797

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn