

RNF185 Rabbit pAb

货号: B18034

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	14kDa/20kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis
细胞定位	Endoplasmic reticulum membrane,Mitochondrion outer membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RNF185 (NP_689480.2).
序列	MASKGPSASASPENSSAGGPSGSSNGAGESGGQDSTFECNICLDTAKDAVISLCGHLFCWPCLHQWLETRPNRQVCPV CKAGISRDKVIPLYGRGSTGQQDPREKTPPRPQGQRPEPENRGGFQGFQFGD

靶点信息

研究背景	E3 ubiquitin-protein ligase that regulates selective mitochondrial autophagy by mediating 'Lys-63'-linked polyubiquitination of BNIP1. Acts in the endoplasmic reticulum (ER-associated degradation (ERAD) pathway, which targets misfolded proteins that accumulate in the endoplasmic reticulum (ER for ubiquitination and subsequent proteasome-mediated degradation. Protects cells from ER stress-induced apoptosis. Responsible for the cotranslational ubiquitination and degradation of CFTR in the ERAD pathway. Preferentially associates with the E2 enzymes UBE2J1 and UBE2J2.
------	--

基因ID	91445
基因名	RNF185
Swiss	Q96GF1
别名	RNF185

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn