

DENND4A Rabbit pAb

货号: B18042

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	209kDa/213kDa
实测分子量	209kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1500-1650 of human DE NND4A (NP_005839.3).
序列	RPGRYFLKSSPSTENMHFPSSISSQTRQSCISTSASGLDTSALSVQGNFDLNSKSKLQENFCTRSIQIPANRSKTAMSKCPI FPMARSISTSGPLDKEDTGRQKLISTGSLPATLQGATDSLGLEWHLPSPDPVTVPYLSPLVVWKELESL

靶点信息

研究背景	This gene encodes a DENN domain-containing protein that may function as a guanine nucleotide exchange factor that specifically activates ras-related protein Rab-10. This protein also contains a interferon stimulated response element-binding domain and may be involved in regulating the v-myc avian myelocytoma tosis viral (MYC) oncogene. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 8.
基因ID	10260

基因名	DENND4A
Swiss	Q7Z401
别名	DENND4A;IRLB;MYCPBP

产品验证



Western blot analysis of DENND4A expressed in HepG2,mouse liver using DENND4A Rabbit pAb at 1:100 0. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn