

CREG1 Rabbit pAb

货号: B18047

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	24kDa
实测分子量	24KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of human CREG1 (NP_003842.1).
序列	LQLSVSNLQENPYATLTMTLAQTNFCKKHGFDQPQSLCVHIMLSGTVTKVNETEMDIAKHSLFIRHPEMKTWPSSHNWFFA

靶点信息

研究背景	The adenovirus E1A protein both activates and represses gene expression to promote cellular proliferation and inhibit differentiation. The protein encoded by this gene antagonizes transcriptional activation and cellular transformation by E1A. This protein shares limited sequence similarity with E1A and binds both the general transcription factor TBP and the tumor suppressor pRb in vitro. This gene may contribute to the transcriptional control of cell growth and differentiation.
基因ID	8804

基因名	CREG1
Swiss	O75629
别名	CREG1;CREG

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn