

RAPGEF1 Rabbit pAb

货号: B18056

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	116kDa/120kDa/122kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Early endosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 844-1095 of human RAP GEF1 (NP_941372.1).
序列	ARGVAARPGLTHDFHSHEIAEQTLTLLDAELFYKIEIPEVLLWAKEQNEEKSPNLTQFTEHFNNMSYWVRSIIMLQEKAQDRE RLLLLKFIKIMKHLRKLNNFNLSYLAISALDSAPIRRLEWQKQTSEGLAEYCTLIDSSSSFRAYRAALSEVEPPCIPYLGILQDLT FVHLGNPDYIDGKVNFSKRWQQFNILDSMRCFQQAHYDMRRNDDIINFNDFSDHLAEEALWELSLKIKPRNITRRKTDR EEKT

靶点信息

研究背景	This gene encodes a human guanine nucleotide exchange factor. It transduces signals from CRK by binding the SH3 domain of CRK, and activating several members of the Ras family of GTPases. This signaling cascade that may be involved in apoptosis, integrin-mediated signal transduction, and cell transformation. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some variants has not been determined.
基因ID	2889
基因名	RAPGEF1
Swiss	Q13905
别名	RAPGEF1;C3G;GRF2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn