

# FREM2 Rabbit pAb

货号: B18076

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                          |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IF/ICC                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                             |
| 理论分子量 | 295kDa/351kDa                                                                                        |
| 实测分子量 | 351kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Jurkat,HepG2,293T                                                                                    |
| 细胞定位  | Cell membrane,Extracellular side,Single-pass type I membrane protein                                 |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                                                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 47-300 of human FREM 2 (NP_997244.4).                                                                                                                                                |
| 序列   | QPAAFGRALLSPGLAGAAGVPAEEAIVLANRGLRVPGFREVWLDPLHDLVLQVQPGDRCAVSVLDNDALAQRPGR LSPK RFPCDFGPGEVRYSHLGARSPSRDRVRLQLRYDAPGGAVVLPLVLEVEVVFTQLEVVTRNLPLVVEELLGTSNALDARSLE FAFQPETEECRVGILSGLGALPRYGELLHYPQVPGGAREGGAPETLLMDCKAFQELGVRYRHTAASRSPNRDWIPMVVELR SRGAPVGSPALKR |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                              |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes an integral membrane protein containing numerous CSPG (chondroitin sulfate proteoglycan element) repeats and Calx-beta domains. The encoded protein localizes to the basement membrane, forming a ternary complex that plays a role in epidermal-dermal interactions. This protein is important for the integrity of skin and renal epithelia. Mutations in this gene are associated with Fraser syndrome. |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |        |
|-------|--------|
| 基因ID  | 341640 |
| 基因名   | FREM2  |
| Swiss | Q5SZK8 |
| 别名    | FREM2  |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)