

LDHD Rabbit pAb

货号: B18088

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	52kDa/54kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,U-87MG,HT-29,Mouse liver
细胞定位	Mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 158-507 of human LDH D (NP_705690.2).
序列	PGADASLCGMAATGASGTNAVRYGTMRDNLNLEVLPDGRLLHTAGRGRHFRFGFWPEIPHHTAWYSPCVSLGRRKS AAGYNLTGLFVGSEGLGLITATTLRLHPAPEATVAATCAFPSVQAQAVDSTVHILQAQVVARIEFLDEVMMDACNRYSKLN CLVAPTLFLEFHGSQQALEEQLRTEEIVQQNGASDFSWAKEAEERSRLWTARHNAWYAALATRPGCKGYSTDVCVPISR LPEIVVQTKEDLNASGLTGSIVGHVGDGNFHCILLVNPDDAEELGRVKAFAEQLGRRALALHGTCTGEHGIGMGKRQLLQ EEVGAVGVETMRQLKAVLDPQGLMNPQKVL

靶点信息

研究背景	The protein encoded by this gene belongs to the D-isomer specific 2-hydroxyacid dehydrogenase family. The similar protein in yeast has both D-lactate and D-glycerate dehydrogenase activities. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.
------	--

基因ID	197257
基因名	LDHD
Swiss	Q86WU2
别名	LDHD;DLD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn