

# AFAP1L2 Rabbit pAb

货号: B18122

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 38kDa/90kDa/91kDa/93kDa                                                                              |
| 实测分子量 | 110kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | U-87MG,LO2                                                                                           |
| 细胞定位  | Cytoplasm                                                                                            |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human AFAP1L 2 (NP_001001936.1).                                                                                                                                                                                                 |
| 序列   | MERYKALEQLLTELDLDFLKILDQENLSSTALVKKSCLAELLRLYTKSSSSDEEYIYMNKVTINKQQNAESQGAPEEQGLLPN<br>GEPSQHSSAPQKSLPDLPPPKMIPERKQLAIPKTESPEGYYEEAEPYDTSLNEDGEAVSSSYESYDEEDGSKGKSAPYQWP<br>SPEAGIELMRDARICAFWLRKKWLGQWAKQLCVIKDNRLLCYKSSKDHSPQLDVNLLGSSVIHKEKQVRKKEHKLKITPM<br>NADVIVLGLQSKDQAEQWLRVIEVSGLPSEGASEGNQYTPDAQRFNCQKPDIAEK |

靶点信息

|      |                                                                                                                                        |
|------|----------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | May play a role in a signaling cascade by enhancing the kinase activity of SRC. Contributes to SRC-regulated transcription activation. |
| 基因ID | 84632                                                                                                                                  |

|       |                                     |
|-------|-------------------------------------|
| 基因名   | AFAP1L2                             |
| Swiss | Q8N4X5                              |
| 别名    | AFAP1L2;CTB-1144G6.4;KIAA1914;XB130 |

产品验证



Western blot analysis of AFAP1L2 expressed in U-87MG,LO2 using AFAP1L2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)