

ZNRF1 Rabbit pAb

货号: B18124

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	23kDa/29kDa
实测分子量	24kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Jurkat,Mouse thymus,Mouse lung,Rat thymus
细胞定位	Cytoplasmic vesicle,Endosome,Lysosome,Membrane,Peripheral membrane protein,secretory vesicle,synaptic vesicle membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human ZNRF1 (NP_115644.1).
序列	MGGKQSTAARSRGPFPGVSTDDSAVPPPGGAPHFGHYRTGGGAMGLRSRSVSSVAGMGM DPSTAGGVPFGLYTPASRGTGDSERAPGGGGASDSTYAHNGYQETGGGHHRDGMLYLGSRASLADALPLHIAPRWFS SHSGFKCPICSKSVASDEMFMHFIMCLSKPRL

靶点信息

研究背景	This gene encodes an E3 ubiquitin-protein ligase that plays a role in neural-cell differentiation. Overexpression of this gene causes neurite-like elongation. The encoded protein contains both a zinc finger and a RING finger motif and is localized in the endosome/lysosome compartment, indicating that it may be involved in ubiquitin-mediated protein modification, and in synaptic vesicle membranes in neurons.
------	--

基因ID	84937
基因名	ZNRF1
Swiss	Q8ND25
别名	ZNRF1;NIN283

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn