

TUBGCP6 Rabbit pAb

货号: B18128

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	196kDa/198kDa/200kDa
实测分子量	250kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse brain
细胞定位	Cytoplasm,centrosome,cytoskeleton,microtubule organizing center
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 600-900 of human TUB GCP6 (NP_065194.2).
序列	KTINLLKLCCPRHYLCWSDVPVPRISVIFSLEELKEIEKDCAVYVGRMERVARHSSVSKEEKLMEIAKQELIAHAREAAASRVLSALSDRQMSERMALDARKREQFQRLKEQFVKDQERRQAARQEELDDDFSYARELRDRERRLKSLEEELERKARQALVDHYSKLSAEAAARREQKALWRIQRHRLSARLRFLLLEDEKHIQEMLKAVSEAHQPQEPDPVLLSVHPQVTSPGPEHPEGQGQCDSGSAEQHSPAWDGWNRPGLLTPQPLKPLAVGAGGRGLQQAEGARPFSDSLSIGDFL

靶点信息

研究背景	The protein encoded by this gene is part of a large multisubunit complex required for microtubule nucleat ion at the centrosome.
基因ID	85378

基因名	TUBGCP6
Swiss	Q96RT7
别名	TUBGCP6;GCP-6;GCP6;MCCRP;MCCRP1;MCPHCR

产品验证



Western blot analysis of TUBGCP6 expressed in A-549,Mouse brain using TUBGCP6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn