

ACTR10 Rabbit pAb

货号: B18164

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	46kDa
实测分子量	46kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,293T,Mouse heart,Mouse kidney,Mouse lung
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 218-417 of human ACTR10 (NP_060947.1).
序列	SDLKRG LKIQA AKFNIDGNNERPSPPPNV DYP LDGEKILHILGSIRDSVVEILFEQDN EEQSVATLILDSLIQCPIDTRKQLAE NLVVIGGTSMLPGFLHRL LAEIRYLVEKPKYKKALG TKTFRIHTPPAKANCVAWLGGAIFGALQDILGSR SVSKEYYNQTGRI PDWCSLNNPPLEMMFDVGKTQPPLMKRAFSTEK

靶点信息

研究背景	Predicted to be involved in retrograde axonal transport of mitochondrion. Predicted to be located in cytosol; extracellular region; and secretory granule. Predicted to be part of dynactin complex.
基因ID	55860

基因名	ACTR10
Swiss	Q9NZ32
别名	ACTR10;ACTR11;Arp10;Arp11;HARP11

产品验证



Western blot analysis of ACTR10 expressed in Mouse heart,SKOV3,293T using ACTR10 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn