

DDX43 Rabbit pAb

货号: B18171

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	21kDa/72kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Mouse thymus,Rat testis
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 429-648 of human DDX 43 (NP_061135.2).
序列	TWPHSVHRLAQSYLKEPMIVYVGTLDLVAVSSVKQNIIVTTEEEKWSHMQTFLQSMSSTDKVIVFVSRKAVADHLSSDLIL GNISVESLHGDREQRDREKALENFKTGKVRILIATDLASRGLDVHDTVHVYNFDFPRNIEEYVHRIGRTGRAGRTGVSITTL TRNDWRVASELINILERANQSIPEELVSMARFKAHQKREMERKMERPQGRPKKFH

靶点信息

研究背景	The protein encoded by this gene is an ATP-dependent RNA helicase in the DEAD-box family and displays tumor-specific expression.
基因ID	55510

基因名	DDX43
Swiss	Q9NXZ2
别名	DDX43;CT13;HAGE

产品验证



Western blot analysis of DDX43 expressed in Mouse testis,Mouse thymus,Rat testis using DDX43 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

实验步骤

访问官网浏览详情: www.ablybio.cn