

RPRD1A Rabbit pAb

货号: B18174

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	31kDa/35kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,HeLa,A-431,Mouse liver,Mouse testis,Rat brain
细胞定位	nucleoplasm,nucleus,RNA polymerase II, holoenzyme
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 47-276 of human RPRD 1A (NP_001290341.1).
序列	PVIVEAFKHVSSETDESCCKHLGRVLSIWEERSVYENDVLEQLKQALYGDKKPRKRTYEQIKVDENENCSSLGSPSEPPQTL DLVRALQDLENAASGDAAVHQRIASLPVEVQEVSLDKITDKESGERLSKMVEDACMLLADYNGRLAAEIDDRKQLTRML ADFLRCQKEALAEKEHKLEEYKRKLARVSLVRKELRSRIQSLPDLRLPNVTGSHMHLPFAGDIYSED

靶点信息

研究背景	This gene encodes a cell-cycle and transcription regulatory protein. The encoded protein interacts with the cell cycle inhibitor cyclin-dependent kinase 4 inhibitor B and may function as a negative regulator of G(1)/S phase progression. This protein also forms homo- and heterodimers with the protein, regulation of nuclear pre-mRNA domain-containing protein 1B, to form a scaffold that interacts with the C-terminal domain of RNA polymerase II subunit B1 and regulates several aspects of transcription. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 16.
------	--

基因ID	55197
基因名	RPRD1A
Swiss	Q96P16
别名	RPRD1A;HsT3101;P15RS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn