

# CDK5RAP1 Rabbit pAb

货号: B18189

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IHC                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:100                                                |
| 理论分子量 | 48kDa/56kDa/57kDa/64kDa/66kDa/67kDa                                                                  |
| 实测分子量 | 68kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | A-549,Mouse thymus                                                                                   |
| 细胞定位  | cytosol,mitochondrion                                                                                |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                                                                               |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human CDK5R AP1 (NP_057492.2).                                                                                                                                                                                                                                                         |
| 序列   | MHPLQCVLQVQRSLGWGPLASVSWLSLRMCRAHSSLSSTMCPSPERQEDGARKDFSSRLAAGPTFQHFLKSASAPQEK<br>LSSEVEDPPPYLMMDELLGRQRKVYLETYGCQMNVDTEIAWSILQKSGYLRTSNLQEADVILLVTC SIREKAEQTIWNRLH<br>QLKALKTRRPRSRVPLRIGILGCMAERLKEEILNREKMVDILAGPDAYRDLPRLLAVAESGQQAANVLLSLDETYADVMPVQ<br>TSASATSAFVSIMRGCDNMCSYCVPFTRGRERSRPIASILEEVKKLSEQGLKEVTLLGQNVNSFRDNSEVQFNSAVPTNLS<br>RGFTTNYKTKQGGLRFAHLLDQVSRV |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                  |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a regulator of cyclin-dependent kinase 5 activity. This protein has also been reported t o modify RNA by adding a methylthio-group and may thus have a dual function as an RNA methylthiotran sferase and as an inhibitor of cyclin-dependent kinase 5 activity. Alternative splicing results in multiple tr anscript variants that encode different isoforms. |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                      |
|-------|--------------------------------------|
| 基因ID  | 51654                                |
| 基因名   | CDK5RAP1                             |
| Swiss | Q96SZ6                               |
| 别名    | CDK5RAP1;C20orf34;C42;CGI-05;HSPC167 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)