

RABGAP1 Rabbit pAb

货号: B18214

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	28kDa/53kDa/67kDa/121kDa
实测分子量	122kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,Mouse brain
细胞定位	Cytoplasm,centrosome,cytoskeleton,cytosol,microtubule organizing center
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 830-1069 of human RAB GAP1 (NP_036329.3).
序列	QQEDPIERFERENRRRLQEANMRLEQENDDLAHELVTSKIALRKDLDNAEEKADALNKELMTKQKLIDAEEEKRRLEEESAQ LKEMCRRELDKAESEIKKNSSIIGDYKQICSQLSERLEKQQTANKVEIEKIRQKVDDCERCREFFNKEGRVKGISSTKEVLDE DTDEEKETLKNQLREMELELAQTKLQLVEAECKIQDLEHHLGLALNEVQAAKKTWFNRTLSSIKTATGVQGKETC

靶点信息

研究背景	May act as a GTPase-activating protein of RAB6A. May play a role in microtubule nucleation by centrosom e. May participate in a RAB6A-mediated pathway involved in the metaphase-anaphase transition.
基因ID	23637

基因名	RABGAP1
Swiss	Q9Y3P9
别名	RABGAP1;GAPCENA;TBC1D11

产品验证



Western blot analysis of RABGAP1 expressed in U-87MG,HeLa,Mouse brain using RABGAP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn