

GCAT Rabbit pAb

货号: B18216

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	45kDa/47kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse eye
细胞定位	Mitochondrion,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 170-419 of human GCA T (NP_055106.1).
序列	KAHKYRYRHLDMADEAKLQEAQKHRLRLVATDGAFSMDGDIAPLQEICCLASRYGALVFMDECHATGFLGPTGRGTDEL LGVMDQVTIINSTLGKALGGASGGYTTGPGPLVSLLRQARPYLFSNSLPPAVVGCASKALDLLMGSNTIVQSMAAKTOR FRSKMEAAGFTISGASHPICPVMLGDARLASRMADDMLKRGIFVIGFSYPVVPKGKARIRVQISAVHSEEDIDRCVEAFVEV GRLHGALP

靶点信息

研究背景	The degradation of L-threonine to glycine consists of a two-step biochemical pathway involving the enzymes L-threonine dehydrogenase and 2-amino-3-ketobutyrate coenzyme A ligase. L-Threonine is first converted into 2-amino-3-ketobutyrate by L-threonine dehydrogenase. This gene encodes the second enzyme in this pathway, which then catalyzes the reaction between 2-amino-3-ketobutyrate and coenzyme A to form glycine and acetyl-CoA. The encoded enzyme is considered a class II pyridoxal-phosphate-dependent aminotransferase. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 14.
基因ID	23464
基因名	GCAT
Swiss	O75600
别名	GCAT;KBL

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn