

METAP1 Rabbit pAb

货号: B18217

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	43kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse kidney,Mouse liver,Mouse pancreas,Rat spleen,Rat testis
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 177-386 of human MET AP1 (NP_055958.2).
序列	RNCYPSPLNYYNFPKSCCTSVNEVICHGIPDRRPLQEGDIVNVDITLYRNGYHGDLNETFFVGEVDDGARKLVQTTYECLM QAIDAVKPGVRYRELGNIIQKHAQANGFSVVRSYCGHGIHKLFTAPNVPHYAKNKAVGVMKSGHVFTIEPMICEGGWQD ETWPDGWTAVTRDGKRSAQFEHTLLVTDTGCEILTRRLDSARPHFMSQF

靶点信息

研究背景	Cotranslationally removes the N-terminal methionine from nascent proteins. The N-terminal methionine i s often cleaved when the second residue in the primary sequence is small and uncharged (Met-Ala-, Cys, Gly, Pro, Ser, Thr, or Val. Required for normal progression through the cell cycle.
基因ID	23173

基因名	METAP1
Swiss	P53582
别名	METAP1;MAP1A;MetAP1A

产品验证



Western blot analysis of METAP1 expressed in Mouse heart,Rat spleen,Mouse kidney using METAP1 Rab bit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30 ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn