

APPBP2 Rabbit pAb

货号: B18232

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	66kDa
实测分子量	67kDa,71kDa,80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,HeLa,Jurkat,Raji,Mouse kidney,Mouse brain,Rat kidney,Rat brain
细胞定位	Cytoplasm,Membrane,Nucleus,Peripheral membrane protein,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 386-585 of human APPB P2 (NP_006371.2).
序列	LEEIAIDCHNKETEQRLLQEADHLHLSSLQLAKKAFGEFNVQTAKHYGNLGRLYQSMRKFKEAEEMHIKAIQIKEQLLGQED YEVALSVGHLASLYNDMNQYENAEKLYLRSIAIGKKLFGEGYSGLEYDYRGLIKLYNSIGNYEKVFHEYHNVLSNWNRLRDR QYSVTDALDVDSTSPQSTEEVVQSFLISQNVGPSC

靶点信息

研究背景	The protein encoded by this gene interacts with microtubules and is functionally associated with beta-am yloid precursor protein transport and/or processing. The beta-amyloid precursor protein is a cell surface p rotein with signal-transducing properties, and it is thought to play a role in the pathogenesis of Alzheimer 's disease. The encoded protein may be involved in regulating cell death. This gene has been found to be highly expressed in breast cancer. Alternative splicing results in multiple transcript variants.
------	---

基因ID	10513
基因名	APPBP2
Swiss	Q92624
别名	APPBP2;APP-BP2;HS.84084;PAT1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn