

CLSTN3 Rabbit pAb

货号: B18244

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	106kDa/107kDa
实测分子量	105kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cell membrane,Endoplasmic reticulum membrane,Golgi apparatus membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 730-850 of human CLSTN3 (NP_055533.2).
序列	SLQQRGLELTNTSAYLTIAGVESITVYEEILRQARYRLRHGAALYTRKFRRLSCSEMNGRYSSENEFIVEVNVLHSMNRVAHPSHVLSSQQFLHRGHQPPPEMAGHSLASSHRNSMIPSAATL

靶点信息

研究背景	May modulate calcium-mediated postsynaptic signals. Complex formation with APBA2 and APP, stabilizes APP metabolism and enhances APBA2-mediated suppression of beta-APP40 secretion, due to the retardation of intracellular APP maturation.
基因ID	9746

基因名	CLSTN3
Swiss	Q9BQT9
别名	CLSTN3;CDHR14;CSTN3;alcbeta

产品验证



Western blot analysis of CLSTN3 expressed in Mouse brain using CLSTN3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn