

FHL5 Rabbit pAb

货号: B18249

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	32kDa
实测分子量	33KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	DU145,PC-3,U-87MG,Mouse heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-284 of human FHL5 (NP_065228.4).
序列	MTTAHFYCQYCTASLLGKKYVLKDDSPYCVTCYDRVFSNYCEECKPIESDSKDLCYKDRHWHEGCFKCTKCNHSLVEKP FAAKDERLLCTECYSNECSSKCFHCKRTIMPGSRKMEFKGNYPWHETCFVCENCRQPIGTKPLISKESGNYCVPCFEKEFAH YCNFCKKVITSGGITFCDQLWHKECFLCSGCRKDLCEEQFMSRDDYPFCVDCYNHLYANKCVACSKPISGLTGAKFICFQD SQWHSECFNCGKCSVSLVGKGLTQNKEIFCQKCGSGMDTDI

靶点信息

研究背景	The protein encoded by this gene is coordinately expressed with activator of cAMP-responsive element m odulator (CREM). It is associated with CREM and confers a powerful transcriptional activation function. CR EM acts as a transcription factor essential for the differentiation of spermatids into mature spermatozoa. There are multiple polyadenylation sites found in this gene. Polymorphisms in this gene may be associate d with susceptibility for migraine headaches. Alternative splicing results in multiple transcript variants en coding the same protein.
------	--

基因ID	9457
基因名	FHL5
Swiss	Q5TD97
别名	FHL5;1700027G07Rik;ACT;FHL-5;dj393D12.2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn