

NSMAF Rabbit pAb

货号: B18258

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	104kDa/107kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NSMAF (NP_003571.2).
序列	MAFIRKKQQEQQLQLYSKERFSLLLLNLEEYFFEQHRANHILHKGSHHERKIRGSLKICSKSVIFEPDSISQPIIKIPLRDCIKIGKHGENGANRHFRTKAKSGGISLIFSQVYFIKEHNVVAPYKIERGKMEYVFELDVPGKVEDVVETLLQLHRASCLDKLGDQTA MITAILQSRLARTSFDKNRFQNISEKLHMECKAEMVTPLVNPGHVCITDTNLYFQPLNGYPKPVVQITLQDVRRIYKRRHGL MPLGLEVFCTEDDLCSDIYLFYEPQDRDDLYFYIATYLEHHVAEHTAES

靶点信息

研究背景	This gene encodes a WD-repeat protein that binds the cytoplasmic sphingomyelinase activation domain of the 55kD tumor necrosis factor receptor. This protein is required for TNF-mediated activation of neutral sphingomyelinase and may play a role in regulating TNF-induced cellular responses such as inflammation . Alternative splicing results in multiple transcript variants.
------	--

基因ID	8439
基因名	NSMAF
Swiss	Q92636
别名	NSMAF;FAN;GRAMD5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn