

TCHH Rabbit pAb

货号: B18266

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	253kDa
实测分子量	280kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Jurakt,HeLa,NIH/3T3,Mouse brain,Mouse kidnay,Mouse lung,Mouse pancreas
细胞定位	cytoskeleton,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human TCHH (NP_009044.2).
序列	RARCDGKESLLQDRRQEEDQRRFEPDRQLEEEPGQRRRQKRQEQERELAEGEEQSEKQERLEQRDRQRRDEELWRQRQEWQEREERRAEEEQLQSCKGHETEEFPDEEQLRRRELLELRKGREEKQQQRRERQDRVFQEEEEKEWRKRETVLRKEEEKLQEEEPQRQRELQEEEEQLRKLERQELRRERQEEEQQQRL

靶点信息

研究背景	The protein encoded by this gene forms crosslinked complexes with itself and keratin intermediate filaments to provide mechanical strength to the hair follicle inner root sheath. The encoded protein also is important for structural integrity of the filiform papillae of the tongue. Defects in this gene are a cause of uncombable hair syndrome.
基因ID	7062

基因名	TCHH
Swiss	Q07283
别名	TCHH;THH;THL;TRHY;UHS3

产品验证



Western blot analysis of TCHH expressed in Mouse brain,U-87MG,Jurakt using TCHH Rabbit pAb at 1:100 0. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn