

GPRIN1 Rabbit pAb

货号: B18339

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	80kDa/102kDa
实测分子量	115kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Hela,HT-29,Mouse brain,Rat brain
细胞定位	Cell membrane,Cell projection,Lipid-anchor,growth cone
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GPRIN1 (NP_443131.2).
序列	MDTAEDPAWLQLLQKDSSPPGPRPTAFFCPQDGLGAGSSAMRDYCPSSQKASPAPPRHTPDQSPGMESRHRSPSGA GEGASCSDGPRGSLACPSPTCFSPQESPSKETLEAHGASISGTPEATTSGKPEPVSSVKTEPKSSDDRNPMFLEKMDFKS SKQADSTSIGKEDPGSSRKADPMFTGKAPEILGKGDPVAPGRMDPMTVRKEDLGLGKVDPLCSSKTYTVSPRKEDPG SLRKVDPVSSDKVDPVFPRKEEPRYSGKEHPVSSEKVAPTSAEKVLDLVLSGKRDPGPGSGKADPM

靶点信息

研究背景	May be involved in neurite outgrowth.
基因ID	114787
基因名	GPRIN1

Swiss	Q7Z2K8
别名	GPRIN1;GRIN1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn