

ZNF598 Rabbit pAb

货号: B18354

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	34kDa/97kDa/98kDa
实测分子量	115kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,A-431,Jurkat,HeLa,LO2
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-420 of human ZNF598 (NP_835461.2).
序列	PRCPPELPFSLFGDLEQHMRQHELFCCRLCLQHLQIFTYERKWYSRKDLARHRMQGDPDDTSHRGHPLCKFCDERYLD NDELLKHLRRDHYFCHFCDSGDAQDYSDYAYLREHFREKHFLECEEGRCSTEQFTHAFRTEIDLKAHRTACHSRRAEAR QNRHIDLQFSYAPRHSRRNEGTVGGEDYEEVDYRSRQGRVARAGTRGAQQSRRGWSRYKREEEDREVA AAVRASVAA QQQEEARRSEDQEEGGRPKKEEAAARGPEDPRGPRRSPRTQGEGPGPKETSTNGPVSQEAFSVTG

靶点信息

研究背景	Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This protein and Grb10-interacting GYF protein 2 have been identified as components of the mammalian 4EHP (m4EHP) complex. The complex is thought to function as a translation repressor in embryonic development.
------	--

基因ID	90850
基因名	ZNF598
Swiss	Q86UK7
别名	ZNF598

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn