

USP44 Rabbit pAb

货号: B18365

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	81kDa
实测分子量	81kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Raji,HT-29,Mouse lung,Mouse pancreas,Mouse liver,Rat stomach,Rat kidney
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-270 of human USP4 4 (NP_115523.2).
序列	RTLSAISKQNYHCTTRSGRFLRSMGTGDDSYFLHDGAQSLLQSEDQLYTALWHRRRILMGKIFRTWFEQSPIGRKKQEEP FQEKIVVKREVKKRRQELEYQVKAELSMPPRKSLRLQGLAQSTIIIVSVQVPAQTPASPAKDKVLSTSENEISQKVSDSSV KRRPIVTP

靶点信息

研究背景	The protein encoded by this gene is a protease that functions as a deubiquitinating enzyme. The encoded protein is thought to help regulate the spindle assembly checkpoint by preventing early anaphase onset. This protein specifically deubiquitinates CDC20, which stabilizes the anaphase promoting complex/cyclos ome.
基因ID	84101

基因名	USP44
Swiss	Q9H0E7
别名	USP44

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn