

MTHFSD Rabbit pAb

货号: B18383

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	41kDa/42kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse brain,Mouse kidney,Rat brain,Rat kidney
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-383 of human MTHFS D (NP_001152849.1).
序列	MEPRAVGVSQDIREQIWGYMESQNLADFPRPVHHRIPNFKGSYLACQNIKDLDFARTQEVKVDPDKPLEGVRLLVLQS KKTLLVPTPRLRTGLFNKITPPPGATKDILRKCATSQGVRNYSVPIGLDSRVLVDLVVVGSAVSEKGWRIGKGEYADLEY AMMVSMGAVSKETPVVTIVHDCQVVDIPEELVEEHDITVDYILTPTRVIATGCKRPKPMGITWFKISLEMMEKIPIRLSLRARE QQAGKDVTLQGEHQHLPEPGCQQTVPVLSVGRPPDTPGPETNSMEAAPGSPPGEGAPLAADVYVGNLPGDARVSDLKR ALRELGSVPLRLTWQGPRRRAFLHYPDSAAAQAVSCLQGLRLGTDTLRVALARQQRDK

靶点信息

研究背景	Enables RNA binding activity. Predicted to be active in cytoplasm.
基因ID	64779

基因名	MTHFSD
Swiss	Q2M296
别名	MTHFSD

产品验证



Western blot analysis of MTHFSD expressed in mouse brain,Rat brain,Mouse kidney using MTHFSD Rab bit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30 ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn