

BIN3 Rabbit pAb

货号: B18405

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	23kDa/29kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,A431,HeLa,Mouse kidney,Mouse spleen,Rat kidney
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human BIN3 (NP_061158.1).
序列	MSWIPFKIGQPKKQIVPKTVERDFEREYGKLQQLQEEQTRRLQKDMKKSTDADLAMSKSAVKISLDLLSNPLCEQDQDLLN MVTALDTAMKRMDAFNQEKVNQIKQTVIEPLKKFGSVFPSLNMAVKRREQALQDYRRLQAKVEKYEEKEKTGPVLAKLHQ AREELRPVREDFEAKNRQLLEEMPRFYGSRLDYFQPSFESLIRAQVYYYSEMHKIFGDLSHQLDQPGHSDEQRENEREAKL SELRALSIVADD

靶点信息

研究背景	The product of this gene is a member of the BAR domain protein family. The encoded protein is comprised solely of a BAR domain which is predicted to form coiled-coil structures and proposed to mediate dimerization, sense and induce membrane curvature, and bind small GTPases. BAR domain proteins have been implicated in endocytosis, intracellular transport, and a diverse set of other processes.
------	---

基因ID	55909
基因名	BIN3
Swiss	Q9NQY0
别名	BIN3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn