

CDCA7L Rabbit pAb

货号: B18408

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	46kDa/48kDa/52kDa
实测分子量	52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Jurkat,Mouse spleen,Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CDCA7 L (NP_061189.2).
序列	MELATRYQIPKEVADIFNAPSDDEEFVGFRRDDVPMETLSSEESCDSFDSLESGKQQDVRFHISKYFTEELRRIFIEDTDSETE DFAGFTQSDLNGKTNPEVMVVEDLSDDGKASLVSEEEEEDEEEDKATPRRSRSRRSSIGLRVAFQFPTKKLANKPDKNSS SEQLFSSARLQNEKKTILERKKDCRQVIQREDSTSESEDDSRDESEQESSDALLKRTMNIKENKAMLAQLLAELNSMPDFFP VRTPTSASRKKTVRRRAFSEGQITRRMNPTRSARPPEKFALENFTVSAAKFAEEFYF

靶点信息

研究背景	Plays a role in transcriptional regulation as a repressor that inhibits monoamine oxidase A (MAOA activity and gene expression by binding to the promoter. Plays an important oncogenic role in mediating the full t ransforming effect of MYC in medulloblastoma cells. Involved in apoptotic signaling pathways; May act do wnstream of P38-kinase and BCL-2, but upstream of CASP3/caspase-3 as well as CCND1/cyclin D1 and E2 F1.
------	---

基因ID	55536
基因名	CDCA7L
Swiss	Q96GN5
别名	CDCA7L;JPO2;R1;RAM2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn