

VILL Rabbit pAb

货号: B18431

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	94kDa/95kDa
实测分子量	96kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse eye,Mouse lung,Rat kidney,Rat eye
细胞定位	actin cytoskeleton,cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 637-856 of human VILL (NP_056957.3).
序列	LDKYDIMLLDTWQEIFLWLGEAASEWKEAWWGQEYLKTHPAGRSPATPIVLVKQGHEPPTFIGWFFTWDPYKWTSHPS HKEVVDGSPAAASTISEITAEVNNLRLSRWPGNGRAGAVALQALKGSQDSENDLVRSPKSAGSRTSSSVSSTSATINGG LRREQLMHQAVEDLPEGVDPARREFYLSDSDFQDIFGKSKEEFYSMATWRQRQEKKQLGFF

靶点信息

研究背景	The protein encoded by this gene belongs to the villin/gelsolin family. It contains 6 gelsolin-like repeats and a headpiece domain. It may play a role in actin-bundling.
基因ID	50853

基因名	VILL
Swiss	O15195
别名	VILL

产品验证



Western blot analysis of VILL expressed in Mouse kidney,Rat kidney,Mouse eye using VILL Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn