

ACOT9 Rabbit pAb

货号: B18449

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	43kDa/46kDa/49kDa/50kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,MCF7,HT-29,Mouse brain,Mouse heart,Mouse kidney,Mouse lung,Rat pancreas
细胞定位	Mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 22-250 of human ACOT 9 (NP_001028755.2).
序列	LTQGPQNPKKQGIFHIHEVRDKLREIVGASTNWRDHSVKAMEERKLLHSFLAKSQDGLPPRRMKDSYIEVLLPLGSEPELRE KYLTVQNTVRFGRILEDLDSLGLVICYMHNKIHSKMSPLSIVTALVDKIDMCKKSLSPQDIKFSGHVSWSVGKTSMEVKM QMFQLHGDEFPCVLDTFVMVARDSENKGPAFVNPLIPESPEEEELFRQGELNKGRRIAFSSTSLLK

靶点信息

研究背景	The protein encoded by this gene is a mitochondrial acyl-CoA thioesterase of unknown function. Two tran script variants encoding different isoforms have been found for this gene.
基因ID	23597

基因名	ACOT9
Swiss	Q9Y305
别名	ACOT9;ACATE2;CGI-16;MT-ACT48;MTACT48

产品验证



Western blot analysis of ACOT9 expressed in Mouse brain,Rat pancreas,U-251MG using ACOT9 Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug p er lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s

实验步骤

访问官网浏览详情: www.ablybio.cn