

COLGALT2 Rabbit pAb

货号: **B18455**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	72kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Mouse brain
细胞定位	Endoplasmic reticulum lumen
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 28-230 of human COLGALT2 (NP_055916.1).
序列	RFVAERDSEDDGEEPVVPESPLQSPTVLVAVLARNAAHTLPHFLGCLERLDYPKSRMAIWAATDHNVDNTEIFREWLKNVQRLYHYVEWRPMDPEPESYPDEIGPKHWPTSRFAHVMKLRQAALRTAREKWSYILFIDVDNFLTNPQTLNLLIAENKTIVAPMLESRGLYSNFWCGITPKGFYKRTPDYVQIREWKRTGCF

靶点信息

研究背景	Beta-galactosyltransferase that transfers beta-galactose to hydroxylysine residues of collagen.
基因ID	23127
基因名	COLGALT2

Swiss	Q8IYK4
别名	COLGALT2;C1orf17;GLT25D2

产品验证



Western blot analysis of COLGALT2 expressed in LO2, Mouse brain using COLGALT2 Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn