

PAN2 Rabbit pAb

货号: B18473

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	135kDa
实测分子量	135kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,NIH/3T3,Mouse thymus
细胞定位	Cytoplasm,Nucleus,P-body
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 700-1000 of human PAN 2 (NP_001120932.1).
序列	DFAQVLKRSICLDQNTQAWCDTCEKYQPTIQTRNIRHLPDILVINCEVNSSKEADFWRMQAEVAFKMAVKKKHGGEISKNK EFALADWKELGSPEGVLVCPSEELKNVWLPFSIRMKMTKNKGLDVCNWTGDEMQWGPARAEEEHGVVYVDLMATVV HILDSRTGGSLVAHIKVGETYHQRKEGVTHQQWYLFNDFLIEPIDKHEAVQFDMNWKVPAILYYVKRNLNSRYNLNIKNP IEASVLLAEASLARKQRKTHTTFTIPLMLNEMPQIGDLVGLDAEFVTLNEEEAELRSDGTKSTI

靶点信息

研究背景	This gene encodes a deadenylase that functions as the catalytic subunit of the polyadenylate binding protein dependent poly(A) nuclease complex. The encoded protein is a magnesium dependent 3' to 5' exoribonuclease that is involved in the degradation of cytoplasmic mRNAs. Alternate splicing results in multiple transcript variants.
------	--

基因ID	9924
基因名	PAN2
Swiss	Q504Q3
别名	PAN2;USP52

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn